



writeWALLS

(dry erase wallcovering)

care and maintenance instructions

Before Writing on Dry Erase Surface for the First Time

Before initial use of the writing surface, the material must be thoroughly cleaned with EXPO® cleaner and a soft cloth, and then towel dried. Follow with clean water and a soft cloth. Use another clean towel to dry. The surface must be rinsed to remove residue left from the cleaning agent.

markers

The use of EXPO® Bullet or Bold markers is recommended. The use of other brands may require more frequent cleaning.

maintenance

- 1. Microfiber cloths are recommended as erasers to remove the dry erase markings.
- 2. Change or clean the cloth often, as they become soiled to avoid leaving excessive marker residue on the surface. Microfiber cloths can be washed.
- 3. Felt erasers may also be used. Change erasers often.
- 4. The use of dirty cloths and/or erasers will spread ink around the surface. As needed, rinse the surface with water and dry with a clean towel.
- 5. For cleaning more marker residue build-up, use EXPO® brand dry erase cleaning solution and a soft towel. Once the surface is clean rinse with water and dry with a soft cloth.
- 6. Projection capable items have a matte/low-glare dry erase surface that requires more frequent cleaning than writing only/high gloss material.
- 7. No markings will be permanent, but it may be necessary to use EXPO® cleaner to bring the surfaces back to "like new" condition.

long term care

- 1. Clean the entire surface on a regular basis to insure the ultimate performance and appearance of the product over the long term.
- 2. Do not use abrasive cleaners/cleaning tools, erasers with hard plastic edges/housing or pointed writing instruments on the dry erase surface.

recommended C eaning & maintenance

Care and maintenance will ensure your dry erase surfaces are clean. Recommended products listed below.

OXIVIR ® TB



One-Step Hospital Grade Disinfectant Cleaner based on proprietary Accelerated Hydrogen Peroxide (AHP ®) technology to deliver fast, effective disinfection and cleaning performance. In a simple ready to use liquid.

- Hospital Grade Disinfectant
- Non-irritant, non-corrosive, non-bleaching, non-rinsing AHP formula
- Virucidal, Bacteriacidal, Fungicidal, Tuberculocidal
- Effective against 27 micro-organisms including MRSA, VRE
- Kills viruses such as Norovirus, COVID-19, H1N1, Influenza A and Hepatitis C in 1 minute

Microban 24



Initially kills 99.9% of germs, including cold and flu viruses for 24 hours. To sanitize hard, non-po-rous, non-food contact surfaces: (1) Let stand 10 seconds, wipe clean with a damp cloth. (2) Apply and allow air to dry without wiping. Reapply if surface is scrubbed. Microban 24 does not provide 24-hour residual virus protection.

- Eliminates odors
- Effective against S. aureus & E. aerogenes bacteria
- Prevents mold & mildew, repeat use every 7 days
- Effective against viruses: Influenza A H1N1, RSV and COVID-19

Isopropyl Alcohol



Alcohol solutions with at least 70 percent alcohol are effective against viruses on hard surfaces. First, clean the surface with water and detergent. Apply the alcohol solution (do not dilute it) and let it sit on the surface for at least 30 seconds to disinfect.

- Powerful general purpose cleaner
- Antiseptic, Antibacterial, Antiviral disinfectant and skin safe.
- Rapidly evaporates upon use
- Non-corrosive leaving zero residue

